

# Cohort Setup for Thermal Assays

*Hannah Chu*

*2018-06-25*

## Contents

1	Model	1
2	Workflow	1
3	Cohort Setup	2
4	Measure	2

## 1 Model

- *Drosophila melanogaster*
- Mated
- 10 cohorts
- 520 individuals
  - 260 male
  - 260 female

## 2 Workflow




Figure 1: Workflow

- One cohort per day
- 4 different treatments per day for all 13 lines
- 50/50 male and female based on cohort, strain, and treatment
- Position, sex, and treatment all randomized

```

library(DiagrammeR)
mermaid("
graph TD
A[Cohort Days]--> B[Slow Ramp]
A --> C[Fast Ramp]
A --> D[Hardening]
A --> E[Static]
B --> F[13 lines]
C --> G[13 lines]
D --> H[13 lines]
E --> I[13 lines]
F --> J[Female]
G --> K[Female]
H --> L[Female]
I --> M[Female]
F --> N[Male]
G --> O[Male]
H --> P[Male]
I --> Q[Male]"
)

```

### 3 Cohort Setup

Keep 3 rearing vials and 1 “fresh” vial for replicates

- 1 old vial (backup)
- 1 middle vial (backup)
- 1 new vial (backup)
- 1 exp. vial
  - dump the adults and wait for newly eclosed adults
  - flip newly eclosed adults into fresh vial consecutive days
  - keep fresh vials for 5 days before using for tests

### 4 Measure

1. Knockdown times
  - based on thermal assay protocol here
2. Survivorship
  - after flies knockdown, move them into 5mL plastic tubes with food
  - monitor every day for individual fly death